

# CDH18 Antibody



PACO08331

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CDH18

**Uniprot**

Q13634

**Synonyms:**

cadherin 18, type 2; CDH18; CDH14; CDH14L; CDH24; EY-CADHERIN ;

**Immunogen:**

Human CDH18.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A