

# CDKN2AIP Antibody



PACO08347

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CDKN2AIP

**Uniprot**

Q9NXV6

**Synonyms:**

CDKN2A interacting protein; CDKN2AIP; CARF; FLJ20036 ;

**Immunogen:**

Human CDKN2AIP.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A