

# CDKN2C Antibody



PACO08349

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CDKN2C

**Uniprot**

P42773

**Synonyms:**

cyclin-dependent kinase inhibitor 2C (p18, inhibits CDK4); CDKN2C; INK4C; p18; p18-INK4C ;

**Immunogen:**

Human CDKN2C.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A