

CEBPA Antibody



PACO08359

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CEBPA

Uniprot

P49715

Synonyms:

CCAAT/enhancer binding protein (C/EBP), alpha; CEBPA; C/EBP-alpha; CEBP ;

Immunogen:

Human CEBPA.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A