

CENPC Antibody



PACO08362

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CENPC

Uniprot

Q03188

Synonyms:

centromere protein C 1; CENPC1; CENP-C; CENPC; MIF2; hcp-4 ;

Immunogen:

Human CENPC1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A