

CENPL Antibody



PACO08365

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CENPL

Uniprot

Q8N0S6

Synonyms:

centromere protein L; CENPL; C1orf155; CENP-L; FLJ31044; FLJ31786; RP3-383J4.1; dJ383J4.3 ;

Immunogen:

Human CENPL.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A