

CEP72 Antibody



PACO08370

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

CEP72

Uniprot

Q9P209

Synonyms:

centrosomal protein 72kDa; CEP72; FLJ10565; KIAA1519; MGC5307 ;

Immunogen:

Human CEP72.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A