

# CEPT1 Antibody



PACO08371

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:****Protein Background:****Gene ID:**

CEPT1

**Uniprot**

Q9Y6K0

**Synonyms:**

choline/ethanolamine phosphotransferase 1; CEPT1; DKFZp313G0615; MGC45223 ;

**Immunogen:**

Human CEPT1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A