

CES3 Antibody



PACO08375

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CES3

Uniprot

Q6UWW8

Synonyms:

carboxylesterase 3; CES3; ES31; FLJ21736 ;

Immunogen:

Human CES3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A