

# CGNL1 Antibody



PACO08387

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

CGNL1

**Uniprot**

Q0VF96

**Synonyms:**

cingulin-like 1; CGNL1; FLJ14957; JACOP; KIAA1749; MGC138254 ;

**Immunogen:**

Human CGNL1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A