

CHAC2 Antibody



PACO08388

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

CHAC2

Uniprot

Q8WUX2

Synonyms:

ChaC, cation transport regulator homolog 2 (E. coli);

Immunogen:

Human CHAC2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A