

# CHD9 Antibody



PACO08400

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CHD9

**Uniprot**

Q3L8U1

**Synonyms:**

chromodomain helicase DNA binding protein 9; CHD9; AD013; CReMM; KISH2; PRIC320 ;

**Immunogen:**

Human CHD9.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A