

# CHMP1B Antibody



PACO08409

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CHMP1B

**Uniprot**

Q7LBR1

**Synonyms:**

chromatin modifying protein 1B; CHMP1B; C10orf2; C18-ORF2; C18orf2; CHMP1.5; Vps46-2 ;

**Immunogen:**

Human CHMP1B.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A