

# CHMP2B Antibody



PACO08411

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CHMP2B

**Uniprot**

Q9UQN3

**Synonyms:**

chromatin modifying protein 2B; CHMP2B; CHMP2.5; DKFZp564O123; DMT1; VPS2-2; VPS2B ;

**Immunogen:**

Human CHMP2B.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A