

# CHMP5 Antibody



PACO08412

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CHMP5

**Uniprot**

Q9NZZ3

**Synonyms:**

chromatin modifying protein 5; CHMP5; C9orf83; CGI-34; HSPC177; PNAS-2; SNF7DC2 ;

**Immunogen:**

Human CHMP5.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A