

CHMP6 Antibody



PACO08413

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CHMP6

Uniprot

Q96FZ7

Synonyms:

chromatin modifying protein 6; CHMP6; FLJ11749; VPS20 ;

Immunogen:

Human CHMP6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A