

CHRFAM7A Antibody



PACO08419

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CHRFAM7A

Uniprot

Q494W8

Synonyms:

CHRNA7 (cholinergic receptor, nicotinic, alpha 7, exons 5-10) and FAM7A (family with sequence similarity 7A, exons A-E) fusion; CHRFAM7A; CHRNA7; CHRNA7-DR1; D-10; MGC120482; MGC120483 ;

Immunogen:

Human CHRFAM7A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A