

# CHST10 Antibody



PACO08425

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CHST10

**Uniprot**

O43529

**Synonyms:**

carbohydrate sulfotransferase 10; CHST10; HNK-1ST; HNK1ST; MGC17148 ;

**Immunogen:**

Human CHST10.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A