

CHST9 Antibody



PACO08430

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CHST9

Uniprot

Q7L1S5

Synonyms:

carbohydrate (N-acetylgalactosamine 4-O) sulfotransferase 9; CHST9; FLJ42328; GALNAC4ST-2 ;

Immunogen:

Human CHST9.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A