

CHSY1 Antibody



PACO08431

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat, Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CHSY1

Uniprot

Q86X52

Synonyms:

chondroitin sulfate synthase 1; CHSY1; CSS1; KIAA0990 ;

Immunogen:

Human CHSY1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A