

# CHURC1 Antibody



PACO08433

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CHURC1

**Uniprot**

Q8WUH1

**Synonyms:**

churchill domain containing 1; CHURC1; C14orf52; FLJ33064; My015; chch ;

**Immunogen:**

Human CHURC1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A