

CIAPIN1 Antibody



PACO08435

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CIAPIN1

Uniprot

Q6FI81

Synonyms:

cytokine induced apoptosis inhibitor 1; CIAPIN1;2810413N20Rik; Anamorsin; DRE2; PRO0915 ;

Immunogen:

Human CIAPIN1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A