

# CKAP2L Antibody



PACO08443

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CKAP2L

**Uniprot**

Q8IYA6

**Synonyms:**

cytoskeleton associated protein 2-like; CKAP2L; FLJ40629; MGC39683 ;

**Immunogen:**

Human CKAP2L.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A