

CKM Antibody



PACO08446

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CKM

Uniprot

P06732

Synonyms:

creatine kinase, muscle; CKM; CKMM; M-CK ;

Immunogen:

Human CKM.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A