

CLDN3 Antibody



PACO08458

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CLDN3

Uniprot

O15551

Synonyms:

claudin 3; CLDN3; C7orf1; CPE-R2; CPETR2; HRVP1; RVP1 ; claudin3

Immunogen:

Human CLDN3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A