

CLEC10A Antibody



PACO08459

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CLEC10A

Uniprot

Q8IUUN9

Synonyms:

C-type lectin domain family 10, member A; CLEC10A; CD301; CLECSF13; CLECSF14; HML; HML2 ;

Immunogen:

Human CLEC10A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A