

CLUAP1 Antibody



PACO08475

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

CLUAP1

Uniprot

Q96AJ1

Synonyms:

clusterin associated protein 1; CLUAP1; FLJ13297; KIAA0643 ;

Immunogen:

Human CLUAP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A