

CNIH3 Antibody



PACO08483

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

CNIH3

Uniprot

Q8TBE1

Synonyms:

cornichon homolog 3 (Drosophila); CNIH3; FLJ38993 ;

Immunogen:

Human CNIH3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A