

CNOT3 Antibody



PACO08493

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CNOT3

Uniprot

O75175

Synonyms:

CCR4-NOT transcription complex, subunit 3; CNOT3; KIAA0691; LENG2; NOT3; NOT3H ;

Immunogen:

Human CNOT3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A