

CNOT8 Antibody



PACO08496

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CNOT8

Uniprot

Q9UFF9

Synonyms:

CCR4-NOT transcription complex, subunit 8; CNOT8; CAF1; CALIF; POP2; hCAF1 ;

Immunogen:

Human CNOT8.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A