

# CNR1 Antibody



PACO08497

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CNR1

**Uniprot**

P21554

**Synonyms:**

cannabinoid receptor 1 (brain); CNR1; CANN6; CB-R; CB1; CB1A; CB1K5; CB1R; CNR ;

**Immunogen:**

Human CNR1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A