

COG8 Antibody



PACO08507

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

COG8

Uniprot

Q96MW5

Synonyms:

component of oligomeric golgi complex 8; COG8; DOR1; FLJ22315 ;

Immunogen:

Human COG8.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A