

COL23A1 Antibody



PACO08509

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

COL23A1

Uniprot

Q86Y22

Synonyms:

collagen, type XXIII, alpha 1; COL23A1; DKFZp434K0621 ;

Immunogen:

Human COL23A1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A