

# COL4A3BP Antibody



PACO08512

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

COL4A3BP

**Uniprot**

Q9Y5P4

**Synonyms:**

collagen, type IV, alpha 3 (Goodpasture antigen) binding protein; COL4A3BP; CERT; CERTL; FLJ20597; GPBP; STARD11 ;

**Immunogen:**

Human COL4A3BP.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A