

# COPS7A Antibody



PACO08528

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

COPS7A

**Uniprot**

Q9UBW8

**Synonyms:**

COP9 constitutive photomorphogenic homolog subunit 7A (Arabidopsis); COPS7A; CSN7A; MGC110877 ;

**Immunogen:**

Human COPS7A.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A