

COPZ1 Antibody



PACO08530

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

COPZ1

Uniprot

P61923

Synonyms:

coatomer protein complex, subunit zeta 1; COPZ1; CGI-120; COPZ; zeta1-COP ;

Immunogen:

Human COPZ1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A