

# COX16 Antibody



PACO08543

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

COX16

**Uniprot**

Q9P0S2

**Synonyms:**

COX16 cytochrome c oxidase assembly homolog (S. cerevisiae); COX16; C14orf112; HSPC203 ;

**Immunogen:**

Human COX16.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A