

COX17 Antibody



PACO08544

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

COX17

Uniprot

Q14061

Synonyms:

COX17 cytochrome c oxidase assembly homolog (S. cerevisiae); COX17; MGC104397; MGC117386 ;

Immunogen:

Human COX17.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A