

COX5B Antibody



PACO08548

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

COX5B

Uniprot

P10606

Synonyms:

cytochrome c oxidase subunit Vb; COX5B; COXVB ;

Immunogen:

Human COX5B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A