

# COX7C Antibody



PACO08556

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

COX7C

**Uniprot**

P15954

**Synonyms:**

cytochrome c oxidase subunit VIIc;

**Immunogen:**

Human COX7C.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A