

COX8A Antibody



PACO08557

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

COX8A

Uniprot

P10176

Synonyms:

cytochrome c oxidase subunit 8A (ubiquitous); COX8A; COX; COX8; COX8-2; COX8L; VIII; VIII-L ;

Immunogen:

Human COX8A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A