

CPA5 Antibody



PACO08561

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CPA5

Uniprot

Q8WXQ8

Synonyms:

carboxypeptidase A5; CPA5; FLJ35969 ;

Immunogen:

Human CPA5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A