

CPLX4 Antibody



PACO08568

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

CPLX4

Uniprot

Q7Z7G2

Synonyms:

complexin 4; CPLX4; CPX-IV; CPXIV; DKFZp686A0185; DKFZp686O0683; FLJ41190; MGC125769; MGC125783; MGC125784 ;

Immunogen:

Human CPLX4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A