

CPNE8 Antibody



PACO08574

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CPNE8

Uniprot

Q86YQ8

Synonyms:

copine VIII; CPNE8; MGC129645; MGC129646 ;

Immunogen:

Human CPNE8.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A