

CPSF1 Antibody



PACO08577

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CPSF1

Uniprot

Q10570

Synonyms:

cleavage and polyadenylation specific factor 1, 160kDa; CPSF1; CPSF160; HSU37012; P/cl.18;

Immunogen:

Human CPSF1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A