

CPSF2 Antibody



PACO08578

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

CPSF2

Uniprot

Q9P2I0

Synonyms:

cleavage and polyadenylation specific factor 2, 100kDa; CPSF2; CPSF100; KIAA1367 ;

Immunogen:

Human CPSF2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A