

CREB3L4 Antibody



PACO08593

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CREB3L4

Uniprot

Q8TEY5

Synonyms:

cAMP responsive element binding protein 3-like 4; CREB3L4; AIBZIP; ATCE1; CREB3; CREB4; JAL; hJAL ;

Immunogen:

Human CREB3L4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A