

CRIP2 Antibody



PACO08600

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CRIP2

Uniprot

P52943

Synonyms:

cysteine-rich protein 2; CRIP2; CRIP; CRP2; ESP1 ;

Immunogen:

Human CRIP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A