

CRYBB3 Antibody



PACO08617

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

CRYBB3

Uniprot

P26998

Synonyms:

crystallin, beta B3; CRYBB3; CATCN2; CRYB3; MGC125772; MGC125773; MGC125774 ;
Beta-B3 crystallin; Beta-crystallin B3

Immunogen:

Human CRYBB3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A