

CRYM Antibody



PACO08619

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CRYM

Uniprot

Q14894

Synonyms:

crystallin, mu; CRYM; DFNA40; THBP ;

Immunogen:

Human CRYM.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A